## Oxford Cem IMPLANT

# Introducing Oxford Cem IMPLANT semi-permanent cement for crowns & bridges on implants



### Good retention but easy to remove? Yes, Oxford Cem IMPLANT can!

Oxford Cem IMPLANT is a unique dual cure luting cement, which combines good fixation with easy removal in case of repair, extra-oral cleaning or peri-implantitis diseases. Oxford Cem IMPLANT is particularly easy to clean up and easy to re-cement.

Features and Benefits	
Features	Benefits
<ul> <li>Controlled slight expansion and low film thickness</li> </ul>	<ul> <li>Safe fixation, excellent marginal seal</li> </ul>
■ Dual cure	<ul> <li>Fast positioning with command set</li> </ul>
<ul> <li>Elastic and shock absorbing properties</li> </ul>	<ul> <li>Reduces bone stress &amp; gives more comfort to patient</li> </ul>
<ul> <li>No bonding, priming or etching</li> </ul>	<ul> <li>Removable clean and easy to re-cement</li> </ul>
<ul> <li>Convenient gel stage</li> </ul>	<ul> <li>Easy removal of excess material</li> </ul>
<ul> <li>Optimised tooth-coloured translucency</li> </ul>	<ul> <li>Esthetic results of high quality</li> </ul>
<ul><li>Eugenol - free</li></ul>	<ul> <li>Compatible with other materials</li> </ul>
<ul> <li>Contains Zincoxide</li> </ul>	<ul> <li>Zinc is known for good antibacterial properties</li> </ul>






### Oxford Cem IMPLANT

#### Oxford Cem IMPLANT: Low stress

Oxford Cem IMPLANT allows a perfect fixation of the restoration on the implant by a small and controlled expansion. At the same time crowns and bridges can be removed easily because of the very low compressive strength and optimal elastic behaviour of the material. Delivered in an automix double syringe, it is a reliable and convenient product.



\*Data available upon request

63

Prl

 $\int$ 

better

### **Ordering Information**

#### Oxford Cem IMPLANT, semi-permanent cement for crowns & bridges on implants

Includes: 5 ml / 8 g Oxford Cem IMPLANT 4:1 I 10 x Oxford Mix TIP(S) for Minimix 4:1 / 10:1, Long

Part Number 12-008 5 ml / 8 g syringe Oxford Cem IMPLANT 4:1 I 10 x Oxford Mix TIP(S) for Minimix 4:1 / 10:1, Long 01-009 50 x Oxford Mix TIP(S) for Minimix 4:1 / 10:1, Long

© Copyright First Scientific Dental Materials GmbH V3 - 01.09.16