


# Oxford GI FILL - Oxford GI CEM

## Oxford GI Fill: Good Mechanical Properties


### Good Reliability!

### Excellent Flexibility!

#### Compressive Strength



#### Flexural Strength



\*Data available upon request

It fits with:

## Oxford GI COAT



### Light Cure Protective Varnish for Conventional Glass Ionomer Cements

Oxford GI COAT is designed to protect the surface of glass ionomer cement (e.g. Oxford GI FILL) fillings from exposure to moisture and saliva during the first hours after placement.

Oxford GI COAT is free of methylmethacrylate.

## Ordering Information

### Oxford GI FILL & Oxford GI CEM, Classical Glass Ionomer Cements for Filling or for Luting

Part Number	Article
51-002A2	15 g powder / 8 ml liquid Oxford GI FILL Handmix Shade A2I dosage spoon
51-002A3	15 g powder / 8 ml liquid Oxford GI FILL Handmix Shade A3 I dosage spoon
51-002A35	15 g powder / 8 ml liquid Oxford GI FILL Handmix Shade A3,5 I dosage spoon
51-001	35 g powder / 20 ml liquid Oxford GI CEM Handmix I dosage spoon

### Oxford GI Coat

Part Number	Article
51-007	5 ml Oxford GI COAT bottle